


TruSystem[®] 7000dV Surgical Table

Redefining Patient Positioning in *da Vinci*[®] Surgery



The Trumpf Medical™ TruSystem[®] 7000dV Surgical Table is designed to work with Integrated Table Motion for the *da Vinci*, *Xi*[™] Surgical System from Intuitive Surgical[®]

Through an exclusive partnership, Intuitive Surgical and Trumpf Medical take operating room integration to a new level.

Integrated Table Motion

Integrated Table Motion enables the *da Vinci® Xi™* Surgical System to connect to Trumpf Medical's TruSystem® 7000dV Operating Table so that a patient can be dynamically positioned while the surgeon operates. The ability to reposition the Operating Table allows for optimal exposure and access to the target anatomy. This helps to increase efficiency and better manage the patient during robotic-assisted procedures.

TruSystem® 7000dV Package

The TruSystem® 7000dV Operating Table is the central piece of a unique package which provides components and accessories for multi-specialty da Vinci Surgery:

- Double Joint Head Section
- Single Part Leg Section
- dV Remote Control
- dV Connection Cable
- Lightweight Transfer Board
- Armboards
- Universal Arm Supports
- Component Cart
- Patient Restraint Strap


With any TruSystem® 7000dV Operating Table two-years of comprehensive service coverage & maintenance are included.



| | | |
|---------------------------|---------------------------|--|
| Patient Weight | 1,000 lbs. (450 kg) | Exceptional weight capacity in any articulation |
| Low Height | 22" ± 0.4" (59 cm ± 1 cm) | Provides ultra low platform for docking access |
| Trendelenburg Tilt | ± 45° ± 30° | Superior flexibility for reaching target anatomy |

Additional Accessories

Trumpf Medical has an extensive range of accessories to compliment any robotic surgery program.



Not all products/ options are available in all countries. For further information about our products or services, please contact your local Trumpf Medical representative or visit our webpage:

www.trumpfmedical.com

Trumpf Medical reserves the right to make changes without notice in design, specifications and models. The only warranty Trumpf Medical makes is the express written warranty extended on the sale or rental of its products. ©2015 TRUMPF Medizin Systeme GmbH + Co. KG, ALL RIGHTS RESERVED.

TRUMPF Medizin Systeme GmbH + Co. KG
Carl-Zeiss-Straße 7-9 · 07318 Saalfeld · Telefon +49 (0) 3671 586-0 · E-Mail info@trumpfmedical.com

February 8th, 2016